Obtaining information from the NMR spectra of complex mixtures: Pseudoseparation and classification by fingerprint


Ricardo Riguera

Center for Research in Biological Chemistry and Molecular Materials (CiQUS)

Universidad de Santiago, Santiago de Compostela, 15782 SPAIN

E-mail: ricardo.riguera@usc.es

The different methods to extract information from the NMR of complex mixtures (industrial, foods, beverages, biological fluids), without previous separation are described. In particular, procedures for selective pseudoseparations by electronic filtering based on relaxation, diffusion or coupling and by the use of selective paramagnetic spin relaxation ions or by combination of those techniques will be presented. Finally, the use of especific mathematical treatment of the data (NMR, MS, GC etc), for the automatic classification of complex mixtures (metabolomic approaches) wil presented and illustrated with examples in the food and medical fields.



References

- 1.- Paramagnetic NMR Relaxation in Polymeric Matrixes: Selective Suppression of Embedded Species and Sensitivity Enhancement. *J. Am. Chem* Soc. **2007**, 129, 15164-15173;
- 2.- NMR methods for unravelling the spectra of complex mixtures. Nat. Prod. Rep., 2011, 28, 78-98;
- 3.- Stepwise Filtering of the Internal Layers of Dendrimers by Transverse Relaxation-Edited NMR. *J. Am. Chem. Soc.* **2013**, **135** (31), 1513-1516
- 4.- The Dynamics of Dendrimers by NMR Relaxation: Interpretation Pitfalls. *J. Am. Chem. Soc.* **2013**, 135, 1972–1977
- 5.- Predicting PSR Filters by Transverse Relaxation Enhancements. Anal. Chem. 2015, 87, 760–767
- 6.- A joint NMR and SPME-GC chemometric approach for very complex mixtures: Grape and zone identification in wines. *Anal. Chem.*, **2016**, 88, 6239–6246.
- 7.- Fast NMR screening of Macromolecular complexes by paramagnetic spin relaxation filter. *ACS Omega* **2018**, *3*, 2974-2983
- 8.-Filtering the NMR Spectra of Complex Mixtures through Polymer-Mediated Paramagnetic Spin Relaxation. *Chem Eur J.* **2018**,24,19236-19242.